

# LEAD BATTERY FOR ASKIR 230/12 BR, 36 BR, C30 BR

Ref. [BCASP165](#)



## Features

CE Marking	CE0123	Medical device	No
Bar code	3661040169095	HS Code	85065090
Warranty (months)	36	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

18/02/2026 21:03:46 – For information purposes only, can be subject to change