


SMALL CHILD CUFF (XS) FOR MOBI

Ref. 513710



SMALL CHILD CUFF (XS) FOR MOBI

Ref. 513710



Description

Sizes available:

- XS arms with a circumference of 9 to 16 cm (513710)
- S arms from 15 to 26 cm in circumference (513711)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446027995	HS Code	90189010
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.143

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1	US Size (XXS - XXL)	XS
------------------------	---	---------------------	----