

## LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513210



CE  
0459



**Light and robust, the blood pressure monitor that combines style and performance.**

### Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

23/02/2026 22:51:23 – For information purposes only, can be subject to change

# **LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)**



**Ref. 513210**

## **Features**

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE0459        | Medical device | No       |
| Bar code          | 3700446009953 | HS Code        | 90189010 |
| Warranty (months) | 24            | Country        | France   |

### **DIMENSIONS AND WEIGHT**

|                           |     |                           |     |
|---------------------------|-----|---------------------------|-----|
| Packed unit / Width (cm)  | 13  | Packed unit / Height (cm) | 7   |
| Packed unit / Length (cm) | 19  | Packed unit / Weight      | 0.5 |
| Bare Unit / Weight (kg)   | 0.4 |                           |     |

### **SPECIFIC STORAGE CONDITIONS**

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 70 |
| Minimum humidity    | 85  | Maximum humidity | 85 |

### **OTHER FEATURES**

|                        |   |  |  |
|------------------------|---|--|--|
| Packed unit / Quantity | 1 |  |  |
|------------------------|---|--|--|