

**LIAN NANO/SILVER RING/ BLACK NYLON /  
MULTI CUFFS (L, M, S)**

**Ref. 513249**



**CE**  
**0459**

GARANTIE  
**2**  
ANS  
2 YEARS  
GUARANTEE

**Light and robust, the blood pressure monitor that combines style and performance.**

## Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

📅 23/02/2026 23:22:59 – For information purposes only, can be subject to change

# LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513249

## Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | CE0459        | Medical device            | No       |
| Bar code                           | 3700446024307 | HS Code                   | 90189010 |
| Warranty (months)                  | 24            | Country                   | France   |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.478    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | -20           | Max temperature           | 70       |
| Minimum humidity                   | 85            | Maximum humidity          | 85       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |

📅 23/02/2026 23:22:59 – For information purposes only, can be subject to change