

**LIAN NANO/SILVER RING/ BLACK NYLON /
MULTI CUFFS (L, M, S)**

Ref. 513260-



CE
0459

GARANTIE
2 ANS
2 YEARS
GUARANTEE

Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

📅 23/02/2026 23:10:45 – For information purposes only, can be subject to change

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513260-

Features

CE Marking	CE0459	Medical device	No
Bar code	3700446007232	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Packed unit / Height (cm)	9
Packed unit / Length (cm)	27	Packed unit / Weight	0.802
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 23/02/2026 23:10:45 – For information purposes only, can be subject to change