

**LIAN CLASSIC / GREY COTTON CUFF /
LARGE ADULT (L)**

Ref. 519001



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE


GARANTIE
3
3 ANS
3 YEARS
GUARANTEE

📅 24/02/2026 23:56:42 – For information purposes only, can be subject to change

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519001



Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

CUFF

- Cotton CUFF, velcro fastening
- Size **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm, ref : 519029), **M** (26-34 cm, ref 519012), or **L** (32-42 cm), **S-M-L** (ref :519050)
- Available in grey, navy, red or patterns
- Supplied with carrying case

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519001

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007225	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.905
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 24/02/2026 23:56:42 – For information purposes only, can be subject to change