

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520622



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 ANS
3 YEARS
GUARANTEE

24/02/2026 23:02:45 – For information purposes only, can be subject to change

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520622



Description

HIGN ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm), **M** (26-34 cm), **L** (32-42 cm) ou **XL** (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520622



Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446006723	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.617
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 24/02/2026 23:02:45 – For information purposes only, can be subject to change