

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

Ref. M1.121.513

LUXAMED®



CE



REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®Ref. M1.121.513

Reflex Hammer according to Buck | small

Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Available in :

White: M1.121.914

Blue :M1.121.812

Pink :M1.121.612

Yellow: M1.121.513

Kid : M1.131.914

Black M1.121.112

Green: M1.121.412

Purple: M1.121.312

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®
Ref. M1.121.513

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960263	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	YELLOW		