

LUER BONE RONGEUR STRAIGHT 16 CM 5 MM

Ref. IP48016



28/02/2026 20:58:35 – For information purposes only, can be subject to change

LUER BONE RONGEUR STRAIGHT 16 CM 5 MM



Ref. IP48016

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040026985	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.119
Bare Unit / Weight (kg)	0.119		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

28/02/2026 20:58:35 – For information purposes only, can be subject to change