


# IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081\_N



02/03/2026 22:24:11 – For information purposes only, can be subject to change

# IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081\_N



Single-head adult stethoscope, designed for reliable and accurate auscultation.

## Description

- Stethoscope with single adult metal horn (Ø: 45 mm)•  
Anti-cold ring and horn of the color of the tube•  
Aluminum lyre incorporated into the tubing
- Also available in double pavilion for adults and a pediatric model

Available in Black (DH7081\_N) or Red (DH7081\_R)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000120	HS Code	9018908400
Warranty (months)	12	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85

### OTHER FEATURES

Colors	Black	Packed unit / Quantity	1
--------	-------	------------------------	---

📅 02/03/2026 22:24:11 – For information purposes only, can be subject to change