

**Ref. GLU010**



**CE**  
**0537**

DISPOSITIF MEDICAL  
**IVD**  
MEDICAL DEVICE

GARANTIE  
**2** ANS  
2 YEARS  
GUARANTEE

**STERILE**



📅 02/03/2026 22:21:49 – For information purposes only, can be subject to change

=====

**Ref. GLU010**



=====

## Description

=====

## Features

CE Marking	CE0537	Medical device	DMDIV (diagnostic In Vitro)
Bar code	4719869705914	HS Code	9027899000
Warranty (months)	24	Sterile	Yes

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Packed unit / Height (cm)	5
Packed unit / Length (cm)	16	Packed unit / Weight	0.335
Bare Unit / Weight (kg)	0.065		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	50
Maximum humidity	90		

### OTHER FEATURES

📅 02/03/2026 22:21:49 – For information purposes only, can be subject to change