

**BOX OF 800 PROBE COVERS FOR TEMPO
EAR PRO AND TEMPO EAR PRO +**

Ref. THE003ACC1



CE



04/03/2026 00:35:22 – For information purposes only, can be subject to change

BOX OF 800 PROBE COVERS FOR TEMPO EAR PRO AND TEMPO EAR PRO +



Ref. THE003ACC1

- Box of 800 probe covers for Tempo EAR PRO and Tempo EAR PRO +

Description

- Box of 800 probe covers for Tempo EAR PRO and Tempo EAR PRO + thermometers
- Probe cover protection for tympanic thermometers
- Dispenser function
- 40 boxes per pack of 20 probe covers
- Single use, for impeccable hygiene
- Designed to be suitable for children and babies
- Color: transparent
- Material : polypropylene

Features

CE Marking	CE	Medical device	Class I
Bar code	6942284605264	HS Code	39269097
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	16
Packed unit / Length (cm)	20	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-25	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES

📅 04/03/2026 00:35:22 – For information purposes only, can be subject to change