

**CA-MI NEW HOSPIVAC 400 SUCTION UNIT,
AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4
VERSION,**



Ref. BCASPI58



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



📅 05/03/2026 22:31:08 – For information purposes only, can be subject to change

CA-MI NEW HOSPIVAC 400 SUCTION UNIT, AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4 VERSION,



Ref. BCASPI58

Description

with foot control

Designed to drain the patient's body fluids, gases, tissues and bones during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm

Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko

plug• Accessories for the FS version:

Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko plug

Pedal

• Accessories for the Full version:

Fluid collection jar with an anti-reflux valve

 05/03/2026 22:31:08 – For information purposes only, can be subject to change

CA-MI NEW HOSPIVAC 400 SUCTION UNIT, AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4 VERSION,



Ref. BCASPI58

- 8x14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12mm
- Antibacterial & hydrophobic filter (single patient)
- Power cord with Schuko plug,
- Pedal
- Electronic Transfer System between collection jars
- Technical characteristics
 - Engine: Oil-free piston pump requiring no maintenance
 - Power supply: 230V-50Hz
 - Maximum vacuum (adjustable): 0.90 Bar -90 kPa -675 mmHg
 - Max free air flow: 90 l/min
 - Noise level: 46.4 dBA
 - Weight: 20 kg
 - Dimensions: 46 x 85 x 42 cm

**CA-MI NEW HOSPIVAC 400 SUCTION UNIT,
AUTOCLAVABLE JARS, 2 X 4000 ML-FS 4
VERSION,**



Ref. BCASPI58

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139449	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	20
Bare Unit / Weight (kg)	20		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 05/03/2026 22:31:08 – For information purposes only, can be subject to change