

# ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99306



## Description

Adult cuff (26 to 35 cm)

Washable nylon cuff without pocket with Velcro fastening

08/03/2026 00:38:33 – For information purposes only, can be subject to change

# ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99306

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040075044	HS Code	9018901000
Warranty (months)	12	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.093	Bare Unit / Weight (kg)	0.093

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 08/03/2026 00:38:33 – For information purposes only, can be subject to change