

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99606



08/03/2026 00:35:11 – For information purposes only, can be subject to change

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99606

Description

Adult cuff (26 to 35 cm)

Washable nylon cuff without pocket with Velcro fastening

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040075310	HS Code	9018901000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.8	Bare Unit / Weight (kg)	0.122
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

08/03/2026 00:35:11 – For information purposes only, can be subject to change