

# SPENGLER ARMBAND, VAQUEZ-LAUBRY NANO POCKET



Ref. 330220

## Description

Made of black Oxford-grade nylon, with inner bladder and dual tubing (10 cm & 60 cm), with split leather manometer holder, and Velcro fastening.

Child/Small adult (S), nylon, velcro, black : 330220

Adult (M), nylon, velcro, black : 330210

Adult (L), nylon, velcro, black : 330200

## Features

CE Marking	NONE	Medical device	No
Bar code	3700446028565	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 09/03/2026 22:05:19 – For information purposes only, can be subject to change