


**POLYESTER AUTOCLAVABLE CUFF  
STRAIGHT WITH 1 BLADDER FOR DDM  
PNEUMATIC TOURNIQUET**

Ref. GBS103

**DDM** 



CE



📅 09/03/2026 22:56:41 – For information purposes only, can be subject to change

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GBS103

Those cuffs are autoclavable!

## Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

Cuff infant (GBS101)	19 x 4 cm
Cuff child (GBS102)	26.5 x 5 cm
Cuff small adult (GBS101)	35.5 x 6 cm
Cuff adult (GBS103)	53 x 6 cm

**POLYESTER AUTOCLAVABLE CUFF  
STRAIGHT WITH 1 BLADDER FOR DDM  
PNEUMATIC TOURNIQUET**



Ref. GBS103

**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203490	HS Code	90189084
Warranty (months)	24	Country	France

**DIMENSIONS AND WEIGHT**

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

**SPECIFIC STORAGE CONDITIONS**

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

**OTHER FEATURES**

09/03/2026 22:56:41 – For information purposes only, can be subject to change