


# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104

**DDM**



CE



09/03/2026 22:13:20 – For information purposes only, can be subject to change

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104



Those cuffs are autoclavable!

## Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

|                           |             |
|---------------------------|-------------|
| Cuff infant (GBS101)      | 19 x 4 cm   |
| Cuff child (GBS102)       | 26.5 x 5 cm |
| Cuff small adult (GBS101) | 35.5 x 6 cm |
| Cuff adult (GBS103)       | 53 x 6 cm   |

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GBS104

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3700468203506 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | France   |

### DIMENSIONS AND WEIGHT

|                           |    |                           |     |
|---------------------------|----|---------------------------|-----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc  |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.5 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES