

**EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)**

**EMERGENCY'S™**

Ref. MEEMER02\_J



GARANTIE  
2 ANS  
2 YEARS  
GUARANTEE



📅 09/03/2026 22:30:25 – For information purposes only, can be subject to change

## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Ref. [MEEMER02\\_J](#)

**Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment**

### Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures

-

- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg

-

Available in waterproof version

- Sold empty

- Available in blue (MEEMER02\_B), yellow (MEEMER02\_J), black (MEEMER02\_N), orange (MEEMER02\_O), pink (MEEMER02\_RC) and red (MEEMER02\_R).

 09/03/2026 22:30:25 – For information purposes only, can be subject to change

# EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)



Ref. [MEEMER02\\_J](#)

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564540782	HS Code	420299000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	26	Bare Unit / Width (cm)	26
Packed unit / Height (cm)	19	Bare Unit / Height (cm)	19
Packed unit / Length (cm)	31	Bare Unit / Length (cm)	31.6
Packed unit / Weight	1.3		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	YELLOW	Packed unit / Quantity	1

📅 09/03/2026 22:30:25 – For information purposes only, can be subject to change