

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE

Ref. 518113



**CE**  
**0459**



11/03/2026 00:28:30 – For information purposes only, can be subject to change

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518113

=====

### Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

# **VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE**



**Ref. 518113**

## **Features**

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                           | 3700446008871 | HS Code                   | 90189010                      |
| Warranty (months)                  | 60            | Country                   | France                        |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.802                         |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | -20           | Max temperature           | 70                            |
| Minimum humidity                   | 20            | Maximum humidity          | 85                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |