

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK

Ref. 507300



12/03/2026 01:10:21 – For information purposes only, can be subject to change

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507300



Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

Availbe in black (507300)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009151	HS Code	9018.90.20.00-9
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	32	Packed unit / Weight	0.393
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

📅 12/03/2026 01:10:21 – For information purposes only, can be subject to change