

# STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507262



CE



12/03/2026 23:47:48 – For information purposes only, can be subject to change

## STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507262



The Magister II is a high-performance professional stethoscope for versatile auscultation.

### Description

#### HAUTE QUALITÉ ACOUSTIQUE

- Membrane flottante multifréquences
- Pression légère = écoute des basses fréquences (vasculaires et cardiaques)
- Pression ferme = écoute des hautes fréquences (pulmonaires)
- Tubulure simple avec un conduit haute performance
- Lyres en aluminium

#### CONFORT

- Bague antifroid
- Embouts auriculaires en silicone
- Nouveau design de lyres plus adapté à l'anatomie
- Indicateur de pavillon

Disponible en version Silver (507251), Black (507261), Myrtille (507262) et Corail (507252)

#### CONDITIONS DE STOCKAGE :

Pression atmosphérique : 700-1060 hPa

Température : -10° à 40° C

Humidité : 95% max

# STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507262



## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042127	HS Code	9018.90.20.00-9
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.425
Bare Unit / Weight (kg)	0.105		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95
<b>OTHER FEATURES</b>			
Colors	Blueberry	Packed unit / Quantity	1

📅 12/03/2026 23:47:48 – For information purposes only, can be subject to change