


## JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. BCASPI77



14/03/2026 23:49:02 – For information purposes only, can be subject to change

# JARS IN POLICARBONATE (MAX 121°C) 400 ML WITH COVER FOR EMIVAC



Ref. [BCASPI77](#)

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040151076	HS Code	4202121900
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	18.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	27
Packed unit / Weight	0.15	Bare Unit / Weight (kg)	0.1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

14/03/2026 23:49:02 – For information purposes only, can be subject to change