

ELITE BAGS DOCTOR, DOCTOR CASE - BLACK

Ref. MEDOCT03_N

ELITE BAGS



GARANTIE
2
ANS
2 YEARS
GUARANTEE



14/03/2026 23:32:54 – For information purposes only, can be subject to change

ELITE BAGS DOCTOR, DOCTOR CASE - BLACK



Ref. MEDOCT03_N

Essential and professional element of the work of the general practitioner, in the context of his home visits or emergencies

Description

- Key closure- 3

extractable compartments

-

Isothermal amp maker with 31 elastics

-

Extractable transparent padded compartment with Velcro dividers for fragile instruments

-

Mini-fastener

-

Briefcase- Shoulder strap and handle

-

Lanyard for positioning on a trolley- Material:

300D polyester high tenacity fabric

-

Weight: 3.35 kg

- Outdoor measures: 40 x 30 x 14 cm / 20.4 L

- Recommended max. capacity: 4 kg

- Sold empty

ELITE BAGS DOCTOR, DOCTOR CASE - BLACK



Ref. [MEDOCT03_N](#)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540010	HS Code	420212
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	40	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	16	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	33	Bare Unit / Length (cm)	40
Packed unit / Weight	3.35	Bare Unit / Weight (kg)	3.35
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 14/03/2026 23:32:54 – For information purposes only, can be subject to change