

LED EXAMINATION LAMP WITH 5 FOOT STAND - LUXAMED ®

Ref. U1.151.014



14/03/2026 23:56:13 – For information purposes only, can be subject to change

LED EXAMINATION LAMP WITH 5 FOOT STAND - LUXAMED ®

Ref. U1.151.014

Optimum lighting for accurate examinations.

Description

The LED examination light (930mm) is designed to provide optimum illumination using energy-efficient LED technology. It features an aluminium illumination head for efficient thermal convection and heat radiation, with a light intensity of approximately 50,000 Lux at a working distance of 30 cm. The lamp uses a 15 V / 11.25 W LED lamp with 4 high-power LEDs, guaranteeing a life of around 50,000 hours.

The colour temperature is approximately 4,200 K. The gooseneck construction with a 90 cm arm allows the desired position to be maintained.

The lamp has a maximum height of 165 cm and is mounted on a stable metal base with castors, powder-coated to prevent tipping.

It has a multifunctional plug (EU/UK/US).

Available in 6 different colors:

White: U1.111.014

White/ice blue: U1.131.014

White/rapeseed yellow: U1.141.014

White/jade green: U1.151.014

White/delicate lilac: U1.161.014

White/magenta : U1.131.014

LED EXAMINATION LAMP WITH 5 FOOT STAND - LUXAMED ®

Ref. U1.151.014

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961406	HS Code	94052190
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	6.28		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Green
--------	-------