

**IDEAL+ DOUBLE HEAD ADULT, BLUE**

Ref. DH7073\_RC



CE



16/03/2026 23:00:35 – For information purposes only, can be subject to change

# IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073\_RC



## Description

- Stethoscope with double metal horn adult (Ø: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073\_B), Black (DH7073\_N), Red (DH7073\_R) or Pink (DH7073\_RC) for the paediatric version (Ø: 35/30 mm)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000106	HS Code	9018908400
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	25	Packed unit / Weight	0.15
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Colors	Pink	Packed unit / Quantity	1

📅 16/03/2026 23:00:35 – For information purposes only, can be subject to change