


# ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201\_V



CE



17/03/2026 22:51:28 – For information purposes only, can be subject to change

# ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201\_V



## Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 20 x 10 x 2 cm

Available colors:

- Azure- Blue

-

Golden

- Grey

- Red

- Green

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067179	HS Code	90184990
Warranty (months)	12	Country	France

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.2
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20
Packed unit / Weight	0.074	Bare Unit / Weight (kg)	0.073

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Colors	Green	Packed unit / Quantity	1
--------	-------	------------------------	---

📅 17/03/2026 22:51:28 – For information purposes only, can be subject to change