

## ASPIRATION PROBE CH 20 (10 PCS)

Ref. BCASPI200

**CA-MI**



### Description

Tracheobronchial suction probes for CA-MI phlegm aspirators  
Probes for use with CA-MI electric Suction Units: Aspiret / Askir 20 / Askir 30 / Askir C30 / Askir 36 BR / New Askir 36 BR Li-ON / New Askir 230/12v BR / New Aspiret

Single Use Probes  
Box of 10 pieces.

 17/03/2026 23:12:43 – For information purposes only, can be subject to change

# ASPIRATION PROBE CH 20 (10 PCS)

Ref. BCASPI200



## Features

CE Marking	NONE	Medical device	No
Bar code	3661040162072	HS Code	90189084
Warranty (months)	12	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85

### OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--