

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M10500

DDM




CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE



📅 18/03/2026 23:42:48 – For information purposes only, can be subject to change

THE SINGLE USE PRESSURE CUFF EASY FUZE DDM

Ref. M10500



The single use pressure cuff made in France !

Description

The pressure infusion cuffs are medical devices designed for daily use in operating theaters, Intensive Care units or emergency transport vehicles. They can be used to manually accelerate infusion of blood and solutions when gravity alone cannot achieve the desired flow rates.

Easyfuse pressure cuffs are single use and compatible with different sizes of blood/solution bags: 500cc, 1000cc and 3000cc.

Optimized monitoring

- **Single use**, no cleaning and no disinfection
- blood/solution bag clearly visible
- pressure indicator with automatic overpressure regulator
- Controlled decompression facilitated by 3 ways valve
- Résistant as made of TPU

All our pressure infusion cuffs are made in France in our Castelculier (Agen) manufacturing site..

**THE SINGLE USE PRESSURE CUFF EASY
FUZE DDM**

Ref. M10500



Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203810	HS Code	98109084
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

18/03/2026 23:42:48 – For information purposes only, can be subject to change