

**VAQUEZ LAUBRY NANO / VELCRO CUFF /  
BLACK / CHILD (S)**

Ref. 512210



**CE**  
**0459**

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**2** ANS  
2 YEARS  
GUARANTEE

## Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflatio bulb. Supplied in a nylon pouch.

19/03/2026 22:54:38 – For information purposes only, can be subject to change

# VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Ref. 512210

## Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008628	HS Code	90189010
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	20	Packed unit / Weight	0.654
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 19/03/2026 22:54:38 – For information purposes only, can be subject to change