

DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.612

LUXAMED®



19/03/2026 22:54:10 – For information purposes only, can be subject to change

DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.612

LUXAMED[®]

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.612

LUXAMED®

Features

| | | | |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381961710 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| Bare Unit / Weight (kg) | 0.04 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -10 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Colors | Green | | |

 19/03/2026 22:54:10 – For information purposes only, can be subject to change