

**Ref. M25000**



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

24/03/2026 22:57:11 – For information purposes only, can be subject to change

=====

**Ref. M25000**

=====

## Description

Pressure cuffs are medical devices designed for use in operating theatres, intensive care units and emergency transport vehicles. They are used to manually accelerate the infusion of blood and solutions when gravity is not sufficient to achieve the desired flow rates.

The GREY FUSE pressure cuff is compatible with all IV and blood bags, and comes in 3 sizes: 500cc (ref: M20500) , 1000cc (ref: M21000) and 3000cc (ref: M25000).

Optimised monitoring and reliability

Reusable and easy to clean/disinfect

Visualisation of the IV bag thanks to a transparent fabric made up of a screen and a polyamide mesh.

Pressure range from 0 to 300mmHg in 10mmHg increments

Pressure gauge accuracy  $\pm 3$  mmHg

Red zone on dial at 300mmHg to prevent overpressure

Pressure relief controlled by a simple tap.

All our cuff are made in our Castelculier factory near Agen.

**Ref. M25000**

## Features

|                                    |               |                           |            |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking                         | CE            | Medical device            | Class I    |
| Bar code                           | 3700468207276 | HS Code                   | 9018908400 |
| Warranty (months)                  | 24            | Country                   | France     |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |            |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | nc         |
| Bare Unit / Weight (kg)            | 0.52          |                           |            |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |            |
| Minimum temperature                | -10           | Max temperature           | 40         |
| Minimum humidity                   | 30            | Maximum humidity          | 40         |
| <b>OTHER FEATURES</b>              |               |                           |            |

24/03/2026 22:57:11 – For information purposes only, can be subject to change