

Ref. BT103S



29/03/2026 19:42:02 – For information purposes only, can be subject to change

=====

Ref. BT103S

=====

Description

=====

Features

| | | | |
|-------------------|---------------|----------------|------------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700468215684 | HS Code | 9018901000 |
| Warranty (months) | 12 | Country | France |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|-----|
| Minimum temperature | 10 | Max temperature | 70 |
| Minimum humidity | 10 | Maximum humidity | 100 |

OTHER FEATURES

 29/03/2026 19:42:02 – For information purposes only, can be subject to change