

Ref. BT105S



29/03/2026 18:58:01 – For information purposes only, can be subject to change

=====

Ref. BT105S

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215707	HS Code	9018901000
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

 29/03/2026 18:58:01 – For information purposes only, can be subject to change