

# CUFF AND LOWER LIMB CUFF WITH 2 BLADDERS FOR DDM TOURNIQUET

Ref. GBS202



Cuffs and lower limb cuff with 2 bladders for DDM pneumatic tourniquets !

## Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **2 removable silicon bladders** for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

Cuff child	2 x (35.5 x 6 cm)
Cuff adult	2 x ( 59 x 6 cm)
Cuff adult	2 x ( 69 x 9 cm)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203513	HS Code	90189084
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

📅 02/04/2026 02:18:51 – For information purposes only, can be subject to change