

# CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS106



Those cuffs are autoclavable!

## Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **1 removable silicon bladder** for easy autoclaving.

Those ones are designed in a **conical form**

Several sizes are available

### Dimensions of the bladder in cm

Lower limb cuff adult (GCCS104)	69 x 9 cm
Lower limb cuff adult L (GCCS105)	81 x 9 cm
Lower limb cuff adult XL (GCCS106)	85,5 x 12 cm

### Dimensions de la poche en cm

Cuissard adulte (GCCS104)	69 x 9 cm
Cuissard adulte L (GCCS105)	81 x 9 cm
Cuissard adulte XL (GCCS106)	85,5 x 12 cm

📅 02/04/2026 02:02:14 – For information purposes only, can be subject to change

# CONICAL LOWER LIMB CUFF WITH 1 BLADDER FOR DDM TOURNIQUET

Ref. GCCS106



## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203575	HS Code	90189084
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

02/04/2026 02:02:14 – For information purposes only, can be subject to change