


# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS103

# DDM



CE



📅 02/04/2026 01:51:41 – For information purposes only, can be subject to change

# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GCS103

Those lower limb cuffs are autoclavable!

## Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

Lower limb cuff infant	35.5 x 6 cm
Lower limb cuff child	53 x 6 cm
Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85.5 x 12 cm

📅 02/04/2026 01:51:41 – For information purposes only, can be subject to change

**AUTOCLAVABLE POLYESTER LOWER LIMB  
CUFF STRAIGHT 1 BLADDER FOR DDM  
PNEUMATIC TOURNIQUET**



**Ref. GCS103**

**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203650	HS Code	90189084
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

**OTHER FEATURES**

02/04/2026 01:51:41 – For information purposes only, can be subject to change