


# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106



02/04/2026 01:58:31 – For information purposes only, can be subject to change

# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106

Those lower limb cuffs are autoclavable!

## Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

Lower limb cuff infant	35.5 x 6 cm
Lower limb cuff child	53 x 6 cm
Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85.5 x 12 cm

# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203681	HS Code	90189084
Warranty (months)	24	Country	France

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

### OTHER FEATURES

📅 02/04/2026 01:58:31 – For information purposes only, can be subject to change