

# PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM

Ref. GMBU04



📅 02/04/2026 00:36:45 – For information purposes only, can be subject to change

# PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM

Ref. GMBU04

Dedicated to emergency situations, the single-use DDM manual pneumatic tourniquet is completely stand-alone, making it a true French innovation!

## Description

Designed to accompany the patient from the accident site to the emergency room, DDM's single-use pneumatic tourniquet is accessible and easy to use, thanks to its safety valve.

This single-use pneumatic tourniquet is self-sufficient in its composition: bulb, cuff and valve. The injured patient leaves and keeps the tourniquet on until surgery.

With a wide cuff and foam reinforcement inside, pressure is distributed more widely, avoiding skin lesions and intense pain.

Without a manometer, thanks to a mechanical valve, this manual pneumatic tourniquet avoids overpressure, so pressure does not exceed 300mmHg.

This tourniquet also features a label directly on the cuff for recording patient data, such as patient name, date and time of tourniquet application.

- Available in :

Adult 23-33cm (GMBU04)

Child (GMBU02)

XL Adult 38-50cm (GMCU05)

# PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM

Ref. [GMBU04](#)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468200024	HS Code	90189084
Country	France		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

02/04/2026 00:36:45 – For information purposes only, can be subject to change