

# ELITE BAGS - MEDICAL LEATHER BRIEFCASE CLASSY'S

Ref. MECLASS01

**ELITE BAGS**



02/04/2026 01:30:43 – For information purposes only, can be subject to change

# ELITE BAGS - MEDICAL LEATHER BRIEFCASE CLASSY'S



Ref. [MECLASS01](#)

**Medical leather bag with traditional design.**

## Description

With compact size to carry essential diagnostic equipment. A professional touch of elegance and distinction.

- Internal pocket for personal items
- With key clasp
- With handle and shoulder strap
- Material: leather
- Outdoor measures: 36 x 21 x 20 cm / 15.12 L
- Weight: 1.15 kg
- Recommended max. capacity: 3 kg
- Sold empty
- Available in brown (MECLASS01) and blue (MECLASS01)

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564540669	HS Code	4202.11.10
Warranty (months)	24	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	24	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	40	Bare Unit / Length (cm)	36
Packed unit / Weight	1.75	Bare Unit / Weight (kg)	1.15

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	30	Maximum humidity	85

### OTHER FEATURES

Colors	Marron	Packed unit / Quantity	1
--------	--------	------------------------	---

📅 02/04/2026 01:30:43 – For information purposes only, can be subject to change