

# CHROME-PLATED LAUBRY STETHOSCOPE / BLACK

Ref. 507300



02/04/2026 00:32:09 – For information purposes only, can be subject to change

# CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507300



## Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

Availbe in black (507300)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009151	HS Code	9018.90.20.00-9
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	32	Packed unit / Weight	0.393
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

📅 02/04/2026 00:32:09 – For information purposes only, can be subject to change