

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK

Ref. 507308



CE



02/04/2026 01:57:04 – For information purposes only, can be subject to change

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507308



Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

Availbe in black (507300)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026523	HS Code	90181990
Warranty (months)	60		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.393

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Red	Packed unit / Quantity	1
--------	-----	------------------------	---

📅 02/04/2026 01:57:04 – For information purposes only, can be subject to change