

Ref. E1.312.012

LUXAMED®



Description

02/04/2026 01:27:23 – For information purposes only, can be subject to change

=====

Ref. E1.312.012

LUXAMED®

Features

| | | | |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381960393 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | 5 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Size | Size 2 | | |

 02/04/2026 01:27:23 – For information purposes only, can be subject to change