

Ref. E1.443.012

**LUXAMED®**



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 ANS  
2 YEARS  
GUARANTEE

📅 02/04/2026 01:37:00 – For information purposes only, can be subject to change

**Ref. E1.443.012**

## Description

McCOY blades for LUXAMED® laryngoscopes with integrated fibre optics are available in several sizes. The end of the optical fibre can be bent further.

	A (mm)	B(mm)
Size 2 (E1.442.012)	115	18
Size 3 (E1.443.012)	135	18
Size 4 (E1.444.012)	155	18

Laryngoscope blade with integrated fibre optics

Removable blade tip facilitates intubation and reduces the risk of injury.

Made from high-quality stainless steel for long life

Easy to clean and sterilise, as there are no dust traps

Maintenance-free

8,000 fibres for maximum life and light intensity

Optimum illumination with the largest possible cross-section (3 x 7mm)

Illuminance of around 14,000 Lux for 3.7V and 2.5V LEDs

Compatible with all handles conforming to ISO 7376 (eco-responsible standard)

Autoclavable up to 134°C (approx. 4,000 times)

Datamatrix with all important information for traceability

=====

Ref. E1.443.012

**LUXAMED®**

## Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960645	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Size	Size 3		

 02/04/2026 01:37:00 – For information purposes only, can be subject to change