

# DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. GBDM103



02/04/2026 02:18:15 – For information purposes only, can be subject to change

# DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. GBDM103



Monobloc straight arm cuff with one bladder and a ribbon

## Description

Monobloc PU décontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and ribbon to secure the Velcro attachment.

They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

		A	B	C
• New born arm cuff	White tube	3.5	29	38
• Child arm cuff	Pink tube	4.5	35	48
• Small adult arm cuff	Parme tube	6	46	61
• Adult arm cuff	Sky-blue tube	8	46	67
• Large adult arm cuff	Turquoise tube	10	47	70

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

02/04/2026 02:18:15 – For information purposes only, can be subject to change

**DECONTAMINABLE ARM CUFF MONOBLOC  
STRAIGHT NEW BORN TOURNIQUET DDM**

Ref. GBDM103



**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468207856	HS Code	90189084
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.5
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

**OTHER FEATURES**

02/04/2026 02:18:15 – For information purposes only, can be subject to change