


**Ref. TS50003S**



**CE**  
**0459**

GARANTIE  
**3**  
3 ANS  
3 YEARS  
GUARANTEE



📅 01/04/2026 23:31:05 – For information purposes only, can be subject to change

=====

**Ref. TS50003S**



=====

## Description

=====

## Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215448	HS Code	90189010
Warranty (months)	36	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100
<b>OTHER FEATURES</b>			

📅 01/04/2026 23:31:05 – For information purposes only, can be subject to change