


Ref. TS50004S



CE
0459

GARANTIE
3
3 YEARS
3 ANS
GUARANTEE



📅 02/04/2026 01:35:44 – For information purposes only, can be subject to change

=====

Ref. TS50004S



=====

Description

=====

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215455	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100
OTHER FEATURES			

📅 02/04/2026 01:35:44 – For information purposes only, can be subject to change