

Ref. C1.416.114

LUXAMED®



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE



02/04/2026 00:43:36 – For information purposes only, can be subject to change

=====

Ref. C1.416.114

LUXAMED®

=====

Description

Available in several colours:

White (C1.416.914)

Light blue (C1.416.934)

Fuchsia (C1.416.974)

Yellow (C1.416.944)

Kids (C1.816.914)

Black (C1.416.114)

Green (C1.416.954)

Purple (C1.416.964)

=====

Ref. C1.416.114

LUXAMED®

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963417	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.17		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		

 02/04/2026 00:43:36 – For information purposes only, can be subject to change