

S/S KIDNEY DISH 20 X 11.5X 4 CM

Ref. AI20124



02/04/2026 23:10:56 – For information purposes only, can be subject to change

S/S KIDNEY DISH 20 X 11.5X 4 CM

Ref. [AI20124](#)



Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant -
- Autoclavable (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040019529 | HS Code | 90189084 |
| Warranty (months) | 12 | Country | Pakistan |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm) | 14 | Bare Unit / Width (cm) | 11.5 |
| Packed unit / Height (cm) | 10 | Bare Unit / Height (cm) | 4 |
| Packed unit / Length (cm) | 23 | Bare Unit / Length (cm) | 21 |
| Packed unit / Weight | 0.133 | Bare Unit / Weight (kg) | 0.133 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 02/04/2026 23:10:56 – For information purposes only, can be subject to change