

# WALL MOUNT STAND FOR NC3 & NC5

COMEN

Ref. BCMONI06

A sturdy wall mount to keep the Comen NC3/NC5 monitor securely in place and easily accessible.

## Description

This wall mount is designed to hold Comen NC3 and NC5 monitors securely. It offers a stable and practical installation, ideal for optimizing space and ensuring easy access in medical environments.

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040219219	HS Code	9018193010
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

03/04/2026 00:26:52 – For information purposes only, can be subject to change